





**DPHL PREPAREDNESS CALL LOG FOR RULE OUT (R/O) ISOLATE  
 (MUST include LIMS Test Requisition Form)-\*see reverse side for R/O testing\***

**BACTERIAL ISOLATE AND SENTINEL TESTS FOR RULE OUT OR CONFIRMATION**

<input type="checkbox"/> <b><i>Bacillus anthracis</i></b> SBA growth/morphology Gram stain Hemolysis Motility @ TTC (tube) Catalase in covered tube under BSC*	<input type="checkbox"/> <b><i>Francisella tularensis</i></b> SBA growth/morphology CHOC growth/morphology MAC growth/morphology Gram stain Urease Catalase in covered tube under BSC* Oxidase Beta-lactamase Satellite/XV required
<input type="checkbox"/> <b><i>Brucella spp.</i></b> SBA morphology/growth MAC morphology/growth Gram stain Catalase in covered tube under BSC* Urease Oxidase Satellite/XV required	<input type="checkbox"/> <b><i>Yersinia pestis</i></b> SBA morphology/growth MAC morphology/growth Gram stain catalase in covered tube under BSC* Oxidase Urease Indole
<input type="checkbox"/> <b><i>Burkholderia mallei/Burkholderia pseudomallei</i></b> SBA growth/morphology/hemolysis SBA @ 42° growth/morphology MAC growth/morphology Gram stain Oxidase Motility (tube) Indole Polymyxin B or Colistin/B. cepacia or TM KB plate for penicillin and amoxicillin/clavulanic acid (Augmentin)	<input type="checkbox"/> <b>Other pertinent information?</b>

\*Catalase test must be performed in the BSC. Use a small glass tube with parafilm wrapped across the opening of the tube.